

The Complete Pressure Relieving Mattress Range from



At a glance mattress selection chart (Please note: higher risk mattresses are perfectly adequate for low risk users also). **RECOMMENDED ✓✓ APPROPRIATE ✓ UNSUITABLE ✗**

	POSITIONING	BARIATRIC	LIGHTWEIGHT	LOW PROFILE	ENVELOPMENT	CUSHION	SIZE RANGE (in)	SIZE RANGE (cm)	MAXIMUM WEIGHT	INSERT WARRANTY
AT RISK	Rotate and flip	✓	✓	✗	✓✓	Viscotech® Standard	74 x 74 x 6	197 x 88 x 15.2	22.0 Stones (140 kg)	4 Years
HIGH RISK	Rotate and flip	✗	✓	✓✓	✓✓	Viscotech® Mattress Overlay	74 x 34 x 6	187.5 x 87.6 x 8.5	17 Stones (178 kg)	4 Years
	Rotate only	✓	✓	✗	✓✓	Viscotech® Plus	77.5 x 34.5 x 6	197 x 88 x 15.2	22 Stones (140 kg)	6 Years
	Rotate only	✓✓	✓	✗	✓✓	Viscotech® Plus HD	77.5 x 34.5 x 6	197 x 88 x 15.2	30 Stones (191 kg)	6 Years
	Rotate only	✓✓	✓	✗	✓✓	Viscotech® Supreme HD	77.5 x 34.5 x 6	197 x 88 x 15.2	30 Stones (191 kg)	8 Years
VERY HIGH RISK	One way only	✓	✓	✗	✓✓	Viscotech® Supreme	77.5 x 34.5 x 6	197 x 88 x 15.2	22 Stones (140 kg)	8 Years
	Rotate and flip	✗	✓	✓✓	✓✓	ViscoPRO® Overlay	76.8 x 34.3 x 2.8	195 x 87 x 7	20 Stones (127 kg)	3 Years
	Rotate only	✓	✓	✗	✓✓	ViscoPRO® Mattress	76.8/82.7 x 34.3 x 6	210 x 87 x 15	20 Stones + (127 kg+)	3 Years
UP TO GRADE 2	One way only	✓✓	✓	✗	✓✓	Integrity® Static Plus	78 x 34.7 x 5.9	198 x 88 x 15.2	39 Stones (247.7 kg)	2 Years
	Rotate and flip	✓✓	✓✓	✓✓	✓✓	Stryker® Sofcare Overlay	74 x 33 x 2	188 x 83 x 5	25 Stones (158.8 kg)	30 Days
	Rotate only	✓✓	✓✓	✗	✓✓	Stryker® Sofcare System	74 x 35 x 5.9	188 x 89 x 15	25 Stones (158.8 kg)	2 Years
	One way only	✓	✓✓	✗	✓✓	Stryker® P100 System	75 x 35 x 7	190.5 x 88.9 x 17.8	21.25 Stones (135 kg)	1 Year
GRADE 4	One way only	✓✓	✓	✗	✓✓	Integrity® Alternating Mattress	77 x 33.9 x 5.9	196 x 86 x 15	28 Stones (178 kg)	2 Years
	One way only	✓✓	✓	✗	✓✓	Integrity® Extra Long Alternating Mattress	86.6 x 33.8 x 5.9	220 x 86 x 15	28 Stones (178 kg)	2 Years
	One way only	✓✓	✓	✗	✓✓	Integrity® Extra Narrow Alternating Mattress	85 x 31.1 x 5.9	216 x 79 x 15	28 Stones (178 kg)	2 Years
	One way only	✓✓	✗	✗	✓✓	Integrity® HD Alternating Mattress	76.8 x 47.25 x 11.8	195 x 120 x 30	20-50 Stones (127-317 kg)	2 Years
	One way only	✓✓	✗	✗	✓✓	Stryker® Comfortgel	80 x 35 x 7	203.2 x 88.9 x 17.8	35.75 Stones (227 kg)	5 Years
	One way only	✓✓	✗	✗	✓✓	Stryker® IsoFlex LAL	80 x 35 x 6	203.2 x 88.9 x 15.2	35.75 Stones (227 kg)	10 Years

The mattress selection advice in this wallchart is based on manufacturers guidelines and should be used in conjunction with individual client needs and appropriate clinical knowledge

Useful background information on pressure ulcers and current guidance can be found at www.epuap.org or www.nice.org.uk

If you require any assistance or further information please contact us on:

01457 890 980
sales@sumed.co.uk www.sumed.co.uk