

The Complete Pressure Relieving Cushion Range from



Intelli-Gel®
cushions




Intelli-Gel® cannot leak or puncture and will continue to perform effectively even if damaged.

Viscotech®
PRESSURE RELIEVING PRODUCTS



High quality foam and viscoelastic foam technology

ACTION®



Range of cushions incorporating Akton® viscoelastic polymer gel.

Viscopro®



High quality viscoelastic "Memory Foam" technology

FlowForm® Ultra 90



Fluid filled module technology with viscoelastic foam casing that can be used any way up or round

stryker®



Uses 3 layer, air channelling, technology to stop the cushion "bottoming out"

At a glance cushion selection chart (Please note: higher risk cushions are perfectly adequate for low risk users also). RECOMMENDED ✓✓ APPROPRIATE ✓ UNSUITABLE ✗

	POSITIONING	ACTIVE USERS	BARIATRIC	LIGHTWEIGHT	LOW PROFILE	ENVELOPMENT	CUSHION	SIZE RANGE (in)	SIZE RANGE (cm)	MAXIMUM WEIGHT	SPECIALS
AT RISK	✗	✗	✗	✓✓	✗	✓	Viscotech® Standard	16x16 to 20x18	40.64x40.64 to 50.8x45.72	17.0 Stones (108 kg)	Yes
	✗	✗	✓✓	✓✓	✗	✓	Viscotech® Standard HD	16x16 to 20x18	40.64x40.64 to 50.8x45.72	28.0 Stones (178 kg)	Yes
	✓✓	✓✓	✗	✓✓	✓✓	✓	Action® Commuter™	16x16 to 20x20	40.64x40.64 to 50.8x50.8	16.0 Stones (102 kg)*	Yes
	✓✓	✓✓	✗	✗	✓✓	✓	Action® Pilot™	16x16 to 20x20	40.64x40.64 to 50.8x50.8	16.0 Stones (102 kg)*	Yes
	✗	✓✓	✗	✓✓	✓✓	✓	Action® Centurian™	16x16 to 20x20	40.64x40.64 to 50.8x50.8	16.0 Stones (102 kg)*	Yes
	✗	✓✓	✓✓	✓	✓✓	✓✓	Inelli-Gel Low Profile	16x16 to 20x18	40.64x40.64 to 50.8x45.72	30.0 Stones (191 kg)	No
	✓✓	✓✓	✓✓	✗	✗	✓✓	Intelli-Gel Contour	16x16 to 20x18	40.64x40.64 to 50.8x45.72	30.0 Stones (191 kg)	No
HIGH RISK	✗	✗	✗	✓✓	✗	✓✓	Viscotech® Plus	16x16 to 20x18	40.64x40.64 to 50.8x45.72	17.0 Stones (108 kg)	Yes
	✗	✗	✓✓	✓✓	✗	✓✓	Viscotech® Plus HD	16x16 to 20x18	40.64x40.64 to 50.8x45.72	28.0 Stones (178 kg)	Yes
	✓✓	✓✓	✗	✓✓	✗	✓✓	Viscotech® Plus Contour	16x16 to 18x20	40.64x40.64 to 45.72x50.8	17.0 Stones (108 kg)	Yes
	✓✓	✓✓	✓✓	✗	✗	✓✓	Action® Xact® Classic	10x10 to 20x20	25.4x25.4 to 50.8x50.8	28.0 Stones (178 kg)	No
	✗	✓✓	✓✓	✗	✓✓	✓	Action® Professional®	16x16 to 20x20	40.64x40.64 to 50.8x50.8	35.0 Stones (222 kg)	Yes
	✗	✓✓	✗	✓	✓✓	✓	Viscopro®	15x15 to 20x18	38.1x38.1 to 50.8x45.72	15.7 Stones (100 kg)	Yes
	✗	✓✓	✓✓	✓	✓✓	✓✓	Intelli-Gel Low Profile	16x16 to 20x18	40.64x40.64 to 50.8x45.72	30.0 Stones (191 kg)	No
	✓✓	✓✓	✓✓	✗	✗	✓✓	Intelli -Gel Contour	16x16 to 20x18	40.64x40.64 to 50.8x45.72	30.0 Stones (191 kg)	No
VERY HIGH RISK	✓✓	✓✓	✓✓	✗	✗	✓✓	Viscotech® Supreme Contour gel	16x16 to 18x20	40.64x40.64 to 45.72x50.8	28.0 Stones (178 kg)	Yes
	✗	✓✓	✗	✗	✗	✓✓	Viscopro®	15x15 to 21x17	38.1x38.1 to 53.34x43.18	23.6 Stones (150 kg)	Yes
	✗	✓✓	✓✓	✓	✓✓	✓✓	Inelli-Gel Low Profile	16x16 to 20x18	40.64x40.64 to 50.8x45.72	30.0 Stones (191 kg)	No
	✓✓	✓✓	✓✓	✗	✗	✓✓	Intelli-Gel Contour	16x16 to 20x18	40.64x40.64 to 50.8x45.72	30.0 Stones (191 kg)	No
THERAPEUTIC 1-2	✗	✗	✗	✗	✗	✓✓	Flowform® Ultra 90	15x16 to 18x20	38.1x40.64 to 45.72x50.8	20.0 Stones (127 kg)	No
	✓✓	✓✓	✗	✗	✓✓	✓✓	Primenest® 2	16x16 to 20x20	40.64x40.64 to 50.8x50.8	30.0 Stones (191 kg)	No
	✗	✓✓	✓✓	✓	✓✓	✓✓	Inelli-Gel Low Profile	16x16 to 20x18	40.64x40.64 to 50.8x45.72	30.0 Stones (191 kg)	No
	✓✓	✓	✓✓	✗	✗	✓✓	Intelli-Gel Contour	16x16 to 20x18	40.64x40.64 to 50.8x45.72	20.0 Stones (191 kg)	No
	✗	✓✓	✓✓	✓✓	✓✓	✓	Stryker® SofCare Cushion	18x19.6 One Size	47x50 One Size	25st. (158.8kg)	No

*subject to size and body surface area

The cushion selection advice in this wallchart is based on manufacturers guidelines and should be used in conjunction with individual client needs and appropriate clinical knowledge

Useful background information on pressure ulcers and current guidance can be found at www.epuap.org or www.nice.org.uk

If you require any assistance or further information please contact us on:

01457 890 980
sales@sumed.co.uk www.sumed.co.uk